



DePuy Synthes Powered Solution for Minimally Invasive Surgery (MIS) of the Foot

The DePuy Synthes specialized
Power Tools solution is designed
to improve the surgical outcome



PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES




Pen-Shaped Power Tools Solution for Minimally Invasive Foot Surgery

New Dedicated Burr Attachment

with high torque and low speed for MIS Foot Surgery¹

- Ergonomically designed for comfortable grip and use in precise applications²
- Speed range : 940 – 3,750 rpm (maximum speed adjustable with speed regulation function on console)
- High Torque (0.22 Nm) prevents lodging or binding of the burrs³
- Irrigation nozzle to help keep the burr cool during use



2 Wedge Burrs

designed to shape bone⁴

5 Shannon Burrs

designed to cut bone⁴

The flutes and lengths of these instruments are specially designed to support you in the challenging surgical work environment.


- Stainless steel
- Sterile-packed
- Single-use

A minimally invasive surgical procedure may shorten operating times, reduce the length of hospital stay and improve cosmesis by reducing the incision length⁵



Electric Pen Drive or Air Pen Drive

- One modular pen for drilling, burring, sawing, screw driving and k-wire setting applications
- Standard console with irrigation
- Hand or foot switch options to suit your preference



PART NUMBER	DESCRIPTION	IMAGE
05.001.128	Drill / Burr Attachment Straight, for Round Shafts \varnothing 2.35 mm / Gear ratio 16:1	
03.228.011S	Wedge Burr 2.9 x 13.2 x 65 mm, Sterile	
03.228.012S	Wedge Burr 4.3 x 13 x 65 mm, Sterile	
03.228.013S	Shannon Burr Short, Helical Flute 2 x 8 x 65 mm, Sterile	
03.228.014S	Shannon Burr Medium, Straight Flute 2 x 11.8 x 65 mm, Sterile	
03.228.015S	Shannon Burr Medium, Helical Flute 2.2 x 12 x 65 mm, Sterile	
03.228.016S	Shannon Burr Long, Straight Flute 2.2 x 20 x 75 mm, Sterile	
03.228.017S	Shannon Burr Long, Straight Flute 3.1 x 20 x 70 mm, Sterile	



For detailed information about the EPD and APD Systems refer to the System Brochure and Instructions for Use / User's Manual
www.jnjmedicaldevices.com

Drill/Burr Attachment



Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland
Tel: +41 61 965 61 11
www.jnjmedicaldevices.com

Manufactured or distributed by:
Synthes USA, LLC
1101 Synthes Avenue
Monument, CO 80132
USA
To order (USA): 800-327-6887
To order (Canada): 855-946-8999

Shannon and Wedge Burrs



Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo - Germany

Distributed in OUS by:
Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland
Tel: +41 61 965 61 11
www.jnjmedicaldevices.com

Distributed in US by:
The Anspach Effort, Inc.
4500 Riverside Drive
Palm Beach Gardens, FL 33410
USA
To order (USA): 800-327-6887
To order (Canada): 855-946-8999

References

1 DePuy Synthes Power Tools. Internal Design Report. 2019. Ref 0000286406

2 DePuy Synthes Power Tools. Internal Design Report. 2019. Ref 0000286406; Functional and Design Requirements. 2020 Ref 0000267363

3 DePuy Synthes Power Tools. Internal Design Reports. 2019. Ref 0000286406, 2020. Ref 0000292909

4 Instruction for Use GA DePuy Synthes 312735. 2019. Ref 0000292029

5 Maffulli N., Easley M: Minimally Invasive Surgery of The Foot and Ankle. ISBN 978-1-84996-416-6 2001: 458-5

